

Stop Cribbing Now With These Easy Steps:

The Miracle Collar[®] must be adjusted correctly to stop your horse from cribbing. Please remember to check the fit frequently to ensure the collar's effectiveness.







Figure 1

Figure 2

Figure 3

- **1** Standing on the near side of your horse, place the wide bottom piece under the jaw with buckles facing you. The **CENTER TAB** should point forward between the jaws.
- Position the shorter NECK STRAP over the neck, through the mane (or over the bridle path) and buckle when it is just snug enough to keep collar in place (Figure 1). You should always be able to slip a flat hand's width between the neck strap and your horse's neck. NOTE: The neck strap does not apply the pressure to stop cribbing and should be adjusted only snug enough to keep the collar in place.
- 3 Run the longer **BROWBAND STRAP** across the forehead just in front of the ears and buckle moderately tight. As when putting anything new on your horse, watch how your horse reacts to the Miracle Collar's adjustment and make further adjustments one hole at a time. (Figure 2)

With your horse's head in the upright position, adjust the browband strap until it applies just enough pressure to stop your horse from cribbing.

5 To make sure the browband strap is adjusted tight enough not to drop over your horse's face, watch your horse eat off the ground. If the browband strap appears loose, continue to tighten it until it can't slide down over your horse's face and eyes. If necessary, punch additional holes in the strap for more adjustment.

Attention: Check the collar on a daily basis as leather becomes softer and more supple from daily use. The browband strap may need to be adjusted occasionally in the head down position described above. The browband strap should always fit tighter than the neck strap.

6 Check the neck strap again and make sure you can slip a flat hand's width between the neck strap and your horse's neck. (Figure 3)

Run straps through metal keepers on the buckles.

Attention: The browband strap applies the pressure to stop cribbing. Be sure it is properly adjusted at all times. The neck strap simply helps hold the Miracle Collar[®] in place.